## Crystal<sup>®</sup> PURE M1

Dedicated to a small batch, *Crystal*<sup>®</sup> Pure M1 is a first plug-and-play isolated aseptic filling unit. Combining AT-Closed Vial<sup>®</sup> technology and last generation isolator from SKAN, it was designed in tight cooperation with users, to offer ergonomic and detail-oriented process.



### Process

Material introduction

- The material is introduced through the left rapid decontamination airlock to the main chamber.
- The Drug Substance is introduced directly to the main chamber via the aseptic liquid connector, AT-Port<sup>™</sup> System, in order to keep the bulk product outside of the isolator and keep it homogenized and cooled, if needed.

#### Filling – Laser re-sealing – Capping

- In the main chamber, the ready-to-fill AT-Closed Vial<sup>®</sup> is first filled through piercing of the stopper with the special needle.
- The puncture trace is then re-sealed by a 1 second laser shot on the stopper surface.
- A snap-fit plastic cap is placed on the vial with a capping tool, protecting the sterility of the stopper.
- Overall batch status is demonstrated live on the integrated screen in the back of the isolator chamber.

#### Material exit

 The filled vials and other materials can be continuously brought outside of the isolator through the mouse house of the right exit LAF.







www.aseptictech.com



# Crystal<sup>®</sup> PURE M1: Specification

## Decontamination with SKANFOG®

Innovative decontamination process SKANFOG<sup>®</sup>, used in *Crystal*<sup>®</sup> Pure M1, is based on the micro-nebulization of hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>).

This technology allows

- achieving of the Total Kill faster, by direct injection of the  $H_2O_2$  in the chamber;
- achieve the aeration down to 0.5 ppm faster;
- short decontamination cycle:
  < 30 min.</li>



Use of the catalytic converter, breaking the  $H_2O_2$ , allows a possibility of air exhaust directly in the room, decreasing installation requirements.

Aseptic filling of liquid parenterals, including

Applications	- Autologous and allogeneic cell therapy, gene therapy;
	- Intermediate products for cell and gene therapy;
	- Immunotherapies;
	- Hospital pharmacy preparations;
	- Individualized or small batch production.
Dimensions	L 3300 x W 982 x H 2272 mm
Weight	ca. 1200 kg
Materials	Body housing: ABS polymer
	Working chamber: Stainless steel AISI 316L (EN 1.4404) surface roughness $\leq$ 0.8
	μm.
	Filling equipment: Mostly stainless steel AISI 316L, PEEK.
Operation pressure	+60 Pa
H <sub>2</sub> O <sub>2</sub> catalyst	Patented SKAN NANOX®
Bulk connection	Drug Substance is located outside of isolator and passed through an AT-Port™
	System
Filling volume	0.1 ml to 50 ml
Filling accuracy	Typically, 1% (over 1 ml, for water-like viscosity product)
Product Path	All single-use
Control system	Embedded control system with 2 touch screen control panels, and 1 screen in the
	main chamber. Batch report generation.
	Two independent power inlets 220-240 V, 16 A
Utilities	Two independent power inlets 220-240 V, 16 A Compressed air of 6 - 10 bar, 10 Nm3/h, ISO 8573-1: 2010 Class 1.3.1.



www.aseptictech.com

