

# Vial Laser Marking Machine *AT-ETCHER*™

## AT-Closed Vial® laser etching

Optimized for all AT-Closed Vial® batch sizes



### ADD NEW INFORMATION JUST BEFORE SHIPPING

### EXTEND EXPIRY DATE THROUGH STABILITY STUDIES

- Identify your vials, laser etch critical information like batch number product name, patient info and expiry date
- Label your frozen vials at any stage of the process

#### Process

- Manually place the ready-to-fill AT-Closed Vial® on its rotative axis
- Let the AT-ETCHER™ accurately laser engrave the epigraph on the collar of the vial

Cycle time	< 20 sec (1-50ml vials)
Dimensions	309x460x700 mm
Utilities	Electricity only
Materials	Compatible with classical sanitizing agents
Validation	CE certificate

© 2023 – Aseptic Technologies S.A. – Brochure AT-ETCHER – V1.3  
Aseptic Technologies S.A. reserves the right to make any changes to the described equipment and characteristics without notice.